

## Product Data Sheet

---

Product Name: 1,3-Phenylenediamine-4-sulfonic Acid

Cat. No.: GD19949

### Chemical Properties

Cas. No. 88-63-1

Formula C<sub>6</sub>H<sub>8</sub>N<sub>2</sub>O<sub>3</sub>S

M.Wt 188.2


Solubility

Storage

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure



**Caution: Product has not been fully validated for medical applications. For research use only.**

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA