

Product Data Sheet

Product Name: 2-Methoxypyridine-3-boronic acid pinacol ester

Cat. No.: GD16232

Chemical Properties

Cas. No. 532391-31-4

Formula C₁₂H₁₈BN₁O₃

M.Wt 235.087


Solubility

Storage

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure



Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA