

Product Data Sheet

Product Name: Fmoc-D-Phe-OH

Cat. No.: GA10900

Chemical Properties

Cas. No. 86123-10-6

Chemical Name (2R)-2-(9H-fluoren-9-ylmethoxycarbonylamino)-3-phenylpropanoic acid

SMILES C1=CC=C(C=C1)CC(C(=O)O)NC(=O)OCC2C3=CC=CC=C3C4=CC=CC=C24


Formula $C_{24}H_{21}NO_4$ M.Wt 387.4

Solubility DMSO : 100 mg/mL (258.11 mM; Storage
Need ultrasonic)

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure



Background

Fmoc-D-Phe-OH

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA