

Product Data Sheet

Product Name: 5,6-Bis[4-(9,9-dimethyl-9,10-dihydroacridine)phenyl]-2,3-dicyano- pyrazine
Cat. No.: GD05076

Chemical Properties

Cas. No. 1883400-34-7

Formula

M.Wt


Solubility

Storage Store at room temperature argon filling

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure



Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA