

## Product Data Sheet

Product Name: DMT-PEG4 Phosphoramidite

Cat. No.: GC81466

### Chemical Properties

Cas. No. 158041-83-9


Formula  $C_{38}H_{53}N_2O_8P$  M.Wt 696.82

Solubility Storage Store at -20°C, stored under nitrogen

General tips For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT, or blue ice upon request.

Structure



### Background

DMT-PEG4 Phosphoramidite is a phosphoramidite that can be used in the synthesis of oligonucleotides.

**Caution: Product has not been fully validated for medical applications. For research use only.**

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA