

## Product Data Sheet

---

Product Name: Boron Trichloride - Methanol Reagent

Cat. No.: GD05309

### Chemical Properties

Cas. No. 36254-91-8

Formula  $\text{CH}_4\text{BCl}_3\text{O}$

M.Wt 149.212


Solubility

Storage  $-20^\circ\text{C}$  argon filling

General tips For obtaining a higher solubility, please warm the tube at  $37^\circ\text{C}$  and shake it in the ultrasonic bath for a while. Stock solution can be stored below  $-20^\circ\text{C}$  for several months.

Shipping Condition Evaluation sample solution: ship with blue ice. All other available size: ship with RT, or blue ice upon request.

Structure



**Caution: Product has not been fully validated for medical applications. For research use only.**

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: [tech@glpbio.com](mailto:tech@glpbio.com)

Address: 10292 Central Ave. #205, Montclair, CA, USA