

Product Data Sheet

Product Name: Desethyl-SID 7969543

Cat. No.: GC79769

Chemical Properties

Cas. No.

Formula $C_{22}H_{20}N_2O_7$

M.Wt 424.4


Solubility

Storage Store at -20°C

General tips For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT, or blue ice upon request.

Structure



Background

Desethyl-SID 7969543 is an isoquinolinone analog of SID 7969543 and a selective SF-1 (NR5A1) inhibitor. Desethyl-SID 7969543 shows no activity against RORA, VP16 and LRH-1, and exhibits no cytotoxicity at concentrations close to its SF-1 IC₅₀ value [1] [2].

References:

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA

Product Data Sheet

- [1]. Madoux F, et al. Potent, selective and cell penetrant inhibitors of SF-1 by functional ultra-high-throughput screening[J]. Molecular pharmacology, 2008, 73(6): 1776-1784.
- [2]. Roth J, et al. Synthesis of small molecule inhibitors of the orphan nuclear receptor steroidogenic factor-1 (NR5A1) based on isoquinolinone scaffolds. Bioorg Med Chem Lett. 2008;18(8):2628-2632.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA