

Product Data Sheet

Product Name: Boc-D-Phenylalaninol

Cat. No.: GA10284

Chemical Properties

Cas. No. 106454-69-7

Chemical Name tert-butyl N-[(2R)-1-hydroxy-3-phenylpropan-2-yl]carbamate

SMILES O=C(OC(C)(C)C)N[C@H](CC1=CC=CC=C1)CO

Formula $C_{14}H_{21}NO_3$

M.Wt 251.3


Solubility

Storage

General tips For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT, or blue ice upon request.

Structure



Background

Boc-D-Phenylalaninol

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA