

Product Data Sheet

Product Name: PhiKan 083 hydrochloride

Cat. No.: GC36896

Chemical Properties

Cas. No. 1050480-30-2

SMILES CNCC1=CC2=C(C=C1)N(CC)C3=C2C=CC=C3.[H]Cl

Formula $C_{16}H_{19}ClN_2$ M.Wt 274.79

Solubility Water: 2 mg/mL (7.28 mM and warming)
Storage Store at -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.


Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA

Product Data Sheet

Structure



Background

PhiKan 083 hydrochloride is a carbazole derivative, which binds to the surface cavity and stabilizes Y220C (a p53 mutant), with a K_d of 167 μ M[1], and a relative binding affinity (K_d) of 150 μ M in Ln229 cells[3]. K_d : 167 μ M (p53-Y220C)[1], 150 μ M (p53Y220C, in Ln229 cells)[3]

PhiKan 083 is a carbazole derivative, which binds to the surface cavity and stabilizes Y220C (a p53 mutant), with a K_d of 167 μ M[1], shows a relative binding affinity (K_d) of 150 μ M for p53Y220C in Ln229 cells[3]. PhiKan 083 slows down its thermal denaturation rate[2]. PhiKan 083 (125 μ M, 48 hours) reduces the cell viability of engineered variants of Ln229 cells[3]. PhiKan 083 (100 μ M) in combination with doxorubicin (1 μ M) enhances the pro-apoptotic activity in all variants of Ln229 cells (p53wt, p53Y220C, p53G245S, p53R282W)[3]. Cell Viability Assay[1] Cell Line: Ln229, Ln229-p53-wt, Ln229-p53-Y220C, Ln229-p53-G245S, Ln229-p53-R282W cells

[1]. Boeckler FM, et al. Targeted rescue of a destabilized mutant of p53 by an in silico screened drug. Proc Natl Acad Sci U S A. 2008 Jul 29;105(30):10360-5. [2]. Rauf SM, et al. Effect of Y220C mutation

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA

Product Data Sheet

on p53 and its rescue mechanism: a computer chemistry approach. Protein J. 2013 Jan;32(1):68-74.
[3]. Paulmurugan R, et al. A protein folding molecular imaging biosensor monitors the effects of drugs that restore mutant p53 structure and its downstream function in glioblastoma cells. Oncotarget. 2018 Apr 20;9(30):21495-21511.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA