

Product Data Sheet

Product Name: TRAP-6 amide

Cat. No.: GC65529

Chemical Properties

Cas. No. 141923-40-2

Formula C₃₄H₅₆N₁₀O₉

M.Wt 748.87


Solubility Water : 100 mg/mL (133.53 mM;
Need ultrasonic)

Storage Store at -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure



Background

TRAP-6 amide is a PAR-1 thrombin receptor agonist peptide[1].

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA

Product Data Sheet

[1]. Sofia Ramström, et al. Platelet Phosphatidylserine Exposure and Procoagulant Activity in Clotting Whole Blood--Different Effects of Collagen, TRAP and Calcium Ionophore A23187. Thromb Haemost. 2003 Jan;89(1):132-41.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA