

Product Data Sheet

Product Name: Polysarcosine50

Cat. No.: GC91876

Chemical Properties

Cas. No.

Chemical Name N-butyl-2-(methyl((methylamino)methyl)amino)acetamide

SMILES CCCCNC(CN(C)C)C(=O)N


Formula $(C_3H_5NO)_n C_6H_{16}N_2$ M.Wt

Solubility Chloroform: Soluble Storage $-20^{\circ}C$

General tips For obtaining a higher solubility, please warm the tube at $37^{\circ}C$ and shake it in the ultrasonic bath for a while. Stock solution can be stored below $-20^{\circ}C$ for several months.

Shipping Condition Evaluation sample solution: ship with blue ice All other available size: ship with RT, or blue ice upon request.

Structure



Background

Polysarcosine50 is an amine-functionalized hydrophilic polymer.¹ It has been used in the synthesis of catechol-functionalized copolypeptoids that can be used to solubilize manganese oxide nanoparticles in aqueous solvents. Surface modification of liposomes with polysarcosine polymers between 10 and 400 units long have been used to increase circulation half-life.²

References:

- [1]. Gebru, H., Wang, X., Li, Z., et al. Brønsted base mediated one-pot synthesis of catechol-ended amphiphilic polysarcosine-b-poly(N-butyl glycine) diblock copolypeptoids. Pure Appl. Chem. 91(3)363-374(2019).
- [2]. Weber, B., Birke, A., Fischer, K., et al. Solution properties of polysarcosine: From absolute and relative molar mass determinations to complement activation. Macromolecules 51(7)2653–2661(2018).

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA

Product Data Sheet

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA