

## Product Data Sheet

Product Name: 2-Chloro-4,6-dimethylphenyl isocyanate

Cat. No.: GD05102

### Chemical Properties

Cas. No. 124421-12-1

Formula C<sub>9</sub>H<sub>8</sub>ClNO

M.Wt 181.619


Solubility

Storage 2-8°C argon filling

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure



**Caution: Product has not been fully validated for medical applications. For research use only.**

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA