

Product Data Sheet

Product Name: Exendin-3/4 (64-86)

Cat. No.: GC60830

Chemical Properties

Cas. No.

Formula $C_{110}H_{169}N_{29}O_{32}$

M.Wt 2409.75

Solubility Storage -20°C, protect from light

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure

EAVRLFIEWLKNGGPSSGAPPPS

Background

Exendin-3/4 (64-86) is a polypeptide from patent CN 106029087 A. The incretin receptor ligands are derived from multiple skin toxicity Shan exosomes -3 skin of SDGTFTSDLSKQM Di EAVRLFIEWLKNGGPSSGAPPPS.

[1]. CN 106029087 A

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA