

Product Data Sheet

Product Name: 5'-Amino-modifier spacer 9 Phosphoramidite

Cat. No.: GC81504

Chemical Properties

Cas. No. 3048248-99-0


Formula $C_{32}H_{46}N_3O_7P$ M.Wt 615.7

Solubility Storage Store at -20°C

General tips For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT, or blue ice upon request.

Structure



Background

5'-Amino-modifier spacer 9 Phosphoramidite is a modified phosphoramidite monomer that can be used for the synthesis of oligonucleotides.

In Vitro, 5'-Amino-modifier spacer 9 Phosphoramidite can be used for solid-phase synthesis of nucleic acids, introducing an amino group at the 5' end for covalent coupling of functional molecules.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA