

Product Data Sheet

Product Name: HPK1-IN-8

Cat. No.: GC63317

Chemical Properties

Cas. No. 1214561-09-7

Formula $C_{19}H_{17}FN_6O_2S$

M.Wt 412.44


Solubility

Storage Store at -20°C

General tips For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT, or blue ice upon request.

Structure



Background

HPK1-IN-8 is an allosteric, inactive conformation-selective inhibitor of full-length HPK1.

[1]. Wang W, et al. Discovery of an Allosteric, Inactive Conformation-Selective Inhibitor of Full-Length HPK1 Utilizing a Kinase Cascade Assay. Biochemistry. 2021 Oct 19;60(41):3114-3124.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA