

## Product Data Sheet

---

Product Name: Magnesium hydrogen phosphate trihydrate

Cat. No.: GB59214

### Chemical Properties

Cas. No. 7782-75-4

Formula  $\text{H}_7\text{MgO}_7\text{P}$

M.Wt 174.33

Solubility

Storage Store at room temperature, keep dry and cool

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure



**Caution: Product has not been fully validated for medical applications. For research use only.**

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA