

Product Name: PROTAC RIPK degrader-2
Cat. No.: GC60307

Cas. No. 1801547-16-9

SMILES CC1=C(SC=N1)C2=CC=C(CNC([C@H]3N(C([C@H](C(C)(C)C)NC(COCCOCCOCCOCOCOC4=CC5=C(C(NC6=CC=C(SC=N7)C7=C6)=CC=N5)C=C4S(C(C)(C)C)(=O)=O)=O)=O)C[C@H](O)C3)=O)C=C2


Formula	C ₅₂ H ₆₅ N ₇ O ₁₁ S ₃	M.Wt	1060.31
---------	---	------	---------

Solubility	Storage
------------	---------

General tips	For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.
--------------	---

Shipping Condition	Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.
--------------------	--

Structure



PROTAC RIPK degrader-2 is a nonpeptidic PROTAC which potently targets serine-threonine kinase RIPK2 and has highly selective for RIPK2 degradation[1].

[1]. Bondeson DP, et al. Catalytic in vivo protein knockdown by small-molecule PROTACs. Nat Chem Biol. 2015 Aug;11(8):611-7.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@qlpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA

Product Data Sheet

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA