

Product Data Sheet

Product Name: rac-TNT-i

Cat. No.: GC79850

Chemical Properties

Cas. No.

Formula $C_{18}H_{25}BrN_2O_2$ M.Wt 381.31

Solubility DMSO: 50 mg/mL (131.13 mM; Need Storage Store at 4°C, protect from light ultrasonic)

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.


Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA

Product Data Sheet

Structure



Background

rac-TNT-i (rac-NPD3064) is the racomer of TNT-i. TNT-i is an inhibitor targeting M-Sec. TNT-i reversibly inhibits the formation of tunnel nanotubes (TNT) induced by M-Sec. TNT-i can reduce the production of wild-type HIV-1 in macrophages and T cells expressing M-Sec. TNT-i has low cytotoxicity towards macrophages and T cells. TNT-i can be used in studies related to HIV-1 infection.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA