

Product Data Sheet

Product Name: (S,R,S)-AHPC-Me-C8-NHS ester

Cat. No.: GC81018

Chemical Properties

Cas. No.


Formula $C_{37}H_{51}N_5O_8S$ M.Wt 725.89

Solubility DMSO: 50 mg/mL (68.88 mM; Need ultrasonic) Storage Store at 4°C, away from moisture and light

General tips For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution: ship with blue ice All other available size: ship with RT, or blue ice upon request.

Structure



Background

The (S,R,S)-AHPC-Me-C8-NHS ester is an E3 ubiquitin ligase VHL ligand, used to recruit the VHL protein. The (S,R,S)-AHPC-Me-C8-NHS ester can be linked to the target protein ligand through a linker to form a PROTAC.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA