

Product Data Sheet

Product Name: Bromo-PEG7-alcohol

Cat. No.: GC67334

Chemical Properties

Cas. No. 86220-35-1


Formula $C_{14}H_{29}BrO_7$ M.Wt 389.28

Solubility Storage Store at -20°C

General tips For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT, or blue ice upon request.

Structure



Background

Bromo-PEG7-alcohol is a polyethylene glycol (PEG) linker consisting of a bromine group and a terminal hydroxyl group. The PEG linker escalates hydrophilicity. The bromide (Br) can serve as a leaving group in nucleophilic substitution reactions (e.g., amino, mercapto). The terminal hydroxyl group can serve as a site for further derivatization reactions, transforming into reactive functional groups such as azide, alkynyl, and reactive esters. Bromo-PEG7-alcohol is used to synthesize PROTACs, coupling molecules, probes, or other compounds with PEG spacer arms^[1].

References:

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA

Product Data Sheet

[1] An S, Fu L. Small-molecule PROTACs: An emerging and promising approach for the development of targeted therapy drugs. EBioMedicine. 2018;36:553-562.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA