

Product Data Sheet

Product Name: Erepdekinra

Cat. No.: GC72187

Chemical Properties

Cas. No. 2641313-47-3

Formula $C_{88}H_{130}N_{22}O_{22}$ M.Wt 1848.11

Solubility DMSO : 100 mg/mL (54.11 mM; Need ultrasonic; Hygroscopic DMSO has a significant impact on the solubility of product, please use newly opened DMSO); H_2O : < 0.1 mg/mL (insoluble)

Storage Store at -20°C, away from moisture and light

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure



Background

Erepdekinra is an interleukin-17A (IL-17A) receptor antagonist.

References:

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA

Product Data Sheet

[1]. WHO Drug Information. International Nonproprietary Names for Pharmaceutical

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA