

Product Data Sheet

Product Name: Sulfo-SMCC activated SureLight™ Allophycocyanin (APC)

Cat. No.: GC52509

Chemical Properties

Cas. No.

Formula

M.Wt 0


Solubility

Storage 2°C to 8°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure



Background

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA

Product Data Sheet

SMCC-Activated SureLight™ - Allophycocyanin

Dye: SMCC-activated SureLight™ Allophycocyanin

Excitation max. λ : 650 nm

Emission max. λ : 660 nm

Uses: Custom antibody labeling, fluorescence-activated cell sorting, flow cytometry, Luminex, and T-cell labeling

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA