

Product Data Sheet

Product Name: Copper(II) trifluoroacetylacetonate

Cat. No.: GD07322

Chemical Properties

Cas. No. 14324-82-4

Formula $C_{10}H_8CuF_6O_4$

M.Wt 369.705

Solubility

Storage Store at room temperature

General tips For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT, or blue ice upon request.

Structure



Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA