

Product Data Sheet

Product Name: Trimethylolpropane tris(2-methyl-1-aziridinepropionate)

Cat. No.: GD22249

Chemical Properties

Cas. No. 64265-57-2

Formula $C_{24}H_{41}N_3O_6$


M.Wt 467.6

Solubility Storage 2-8°C, Sealed in dry.

General tips For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT, or blue ice upon request.

Structure



Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA