

Product Data Sheet

Product Name: WAY-600868

Cat. No.: GD23683

Chemical Properties

Cas. No. 348-45-8

Formula

M.Wt

Solubility DMSO : 125 mg/mL (511.69 mM;
Need ultrasonic)

Storage -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure



Background

Protein kinase inhibitor 6 is a protein kinase inhibitor.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA