

Product Data Sheet

Product Name: immunoglobulin light chain variable region fragment [Homo sapiens]

Cat. No.: GP10100

Chemical Properties

Cas. No.

SMILES NC(C(NC(C(C)O)C(NC(CC(C)C)C(NC(CCCCN)C(NC(C(C)CC)C(NC(CO)C(NC(CCCNC(N)=N)C(O)=O)=O)=O)=O)=O)=O)CC1=CC=CC=C1


Formula $C_{40}H_{69}N_{11}O_{10}$ M.Wt 864.04

Solubility $\geq 86.4\text{mg/mL}$ in DMSO Storage Store at -20°C

General tips For obtaining a higher solubility, please warm the tube at 37°C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution: ship with blue ice. All other available size: ship with RT, or blue ice upon request.

Structure



Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA