

Product Data Sheet

Product Name: Anti-pan MHC Class II (HLA II) Antibody (IVA12)

Cat. No.: GC79402

Chemical Properties

Cas. No.

Formula

M.Wt

Solubility

Storage Store at -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure



Background

Anti-pan MHC Class II (HLA II) Antibody (IVA12) recognizes human HLA-II (HLA-DR, HLA-DP, and HLA-DQ). Anti-pan MHC Class II (HLA II) Antibody (IVA12) blocks MHC-II in vitro. Recommend

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA

Product Data Sheet

Isotype Controls: Mouse IgG1 kappa, Isotype Control [1].

References:

[1]. Miranda M, et al. FVIII peptides presented on HLA-DP and identification of an A3 domain peptide binding with high affinity to the commonly expressed HLA-DP4. Haematologica. 2025 Jun 1;110(6):1316-1327.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA