

Product Data Sheet

Product Name: Preangiotensinogen (11-14) (human)

Cat. No.: GA23375

Chemical Properties

Cas. No. 102813-98-9

Formula $C_{21}H_{35}N_7O_6$ M.Wt 481.55

Solubility Soluble in DMSO Storage Store at -20°C

General tips For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT, or blue ice upon request.

Structure



Background

Cleavage product of human preangiotensinogen (1-14) (M-1120) by the action of renin.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA