

Product Data Sheet

Product Name: PKC β pseudosubstrate

Cat. No.: GC15471

Chemical Properties

Cas. No. 172308-76-8

[illegible]


Formula	C ₁₇₇ H ₂₉₄ N ₆₂ O ₃₈ S ₃	M.Wt	3994.84
---------	--	------	---------

Solubility	Soluble to 1 mg/ml in Water	Storage	Store at -20°C
------------	-----------------------------	---------	----------------

General tips	For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.
--------------	---

Shipping Condition	Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.
--------------------	--

Structure



Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA