

## Product Data Sheet

Product Name: Fmoc-Threoninol(tBu) DHP HM Resin

Cat. No.: GA10521

### Chemical Properties

Cas. No. N/A

Formula 100-200:mesh,1%DVB

M.Wt 0.3~0.8mmol/g

Solubility

Storage

**General tips** For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

**Shipping Condition** Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure



### Background

Fmoc-Threoninol(tBu) DHP HM Resin

**Caution: Product has not been fully validated for medical applications. For research use only.**

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA