

Product Data Sheet

Product Name: S6 Kinase Substrate Peptide 32

Cat. No.: GP10104

Chemical Properties

Cas. No.


Formula $C_{149}H_{270}N_{56}O_{49}$ M.Wt 3630.1

Solubility $\geq 36.3\text{mg/mL}$ in DMSO Storage Store at -20°C

General tips For obtaining a higher solubility, please warm the tube at 37°C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution: ship with blue ice All other available size: ship with RT, or blue ice upon request.

Structure



Background

Peptide that measures the activity of kinases that phosphorylate ribosomal protein S6. It can also be used as a substrate.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA