

Product Data Sheet

Product Name: 3-Methoxyphenylboronic acid (contains varying amounts of Anhydride)

Cat. No.: GD16079

Chemical Properties

Cas. No. 10365-98-7


Formula C₇H₉BO₃ M.Wt 151.956

Solubility DMSO : 100 mg/mL (658.07 mM; Need ultrasonic) Storage 2-8°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure



Background

3-Methoxybenzeneboronic acid is a biochemical reagent that can be used as a biological material or organic compound for life science related research.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA