

Product Data Sheet

Product Name: Jedi2
Cat. No.: GC50519

Chemical Properties

Cas. No. 651005-90-2

SMILES CC1=C(C(O)=O)C=C(C2=CC=CS2)O1


Formula C₁₀H₈O₃S M.Wt 208.23

Solubility DMSO : 50 mg/mL (240.12 mM; Need ultrasonic) Storage Store at -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure



Background

Piezo1 channel activator. Promotes Ca²⁺ influx into HEK293 cells expressing Piezo1 channels. Elicits responses in mouse Piezo1-transfected, but not Piezo2-transfected cells.

Wang et al (2018) A lever-like transduction pathway for long-distance chemical- and mechano-gating of the mechanosensitive Piezo1 channel. Nat.Commun. 9 1300 PMID:29610524

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA