

Product Data Sheet

Product Name: TNKS-IN-1

Cat. No.: GC79063

Chemical Properties

Cas. No. 59455-93-5

Formula $C_{18}H_{18}N_2O$ M.Wt 278.36

Solubility DMSO: 10 mg/mL (35.92 mM; Need ultrasonic) Storage Store at -20°C

General tips For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT, or blue ice upon request.

Structure



Background

TNKS-IN-1 is a selective and potent tankyrase (TNKS1 / TNKS2) inhibitor with IC₅₀ values of 7.9 nM and 8.7 nM. TNKS-IN-1 induces axin2 accumulation with an EC 50 of 1500 nM. TNKS-IN-1 can antagonize the Wnt signal transduction pathway by stabilizing axin proteins and promoting β-catenin degradation. TNKS-IN-1 can be used for the research of colorectal cancer [1].

In Vitro, TNKS-IN-1 (Compound 3a) (1 h) potently inhibits purified human TNKS1 and TNKS2 with

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA

Product Data Sheet

IC₅₀ values of 7.9 nM and 8.7 nM, respectively, and exhibits >10000 nM for PARP1 [1]. TNKS-IN-1 (24 h) induces axin2 accumulation in DLD1 cells with an EC₅₀ of 1500 nM, confirming cellular tankyrase inhibition activity^[1].

References:

[1]. Buchstaller HP, et al. Discovery and Optimization of 2-Arylquinazolin-4-ones into a Potent and Selective Tankyrase Inhibitor Modulating Wnt Pathway Activity. J Med Chem. 2019;62(17):7897-7909.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA