

Product Data Sheet

Product Name: PH14
Cat. No.: GC73268

Chemical Properties

Cas. No.


Formula $C_{34}H_{39}N_9O_5$ M.Wt 653.73

Solubility DMSO : 50 mg/mL (76.48 mM; Need ultrasonic) Storage -20°C

General tips For obtaining a higher solubility , please warm the tube at 37°C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure



Background

PH14 is a dual PI3K/HDAC inhibitor with IC₅₀ values of 20.3 nM and 24.5 nM for PI3K α and HDAC3, respectively. PH14 has antiproliferative activity and also induces apoptosis in Jeko-1 cells. PH14 can be used in cancer research, such as lymphoma.

References:

[1]. Deng J, et al. Discovery of benzamide-based PI3K/HDAC dual inhibitors with marked pro-apoptosis activity in lymphoma cells. Eur J Med Chem. 2023 Nov 3;262:115915.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA

Product Data Sheet

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA