

Product Data Sheet

Product Name: PAMAM-OH dendrimer, generation 3 solution

Cat. No.: GF04945

Chemical Properties

Cas. No.

Formula $[\text{NH}_2(\text{CH}_2)_2\text{NH}_2]:(\text{G}=3); \text{dendri PAM}$ M.Wt 6941
 $\text{AM}(\text{NHCH}_2\text{CH}_2\text{OH})_{32}$

Solubility Storage 2-8°C

General tips For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT, or blue ice upon request.

Structure



Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA