

Product Data Sheet

Product Name: Ald-Ph-amido-PEG2-C2-Pfp ester

Cat. No.: GC38337

Chemical Properties

Cas. No. 2101206-60-2

SMILES O=C(OC1=C(F)C(F)=C(F)C(F)=C1F)CCOCCOCCNC(C2=CC=C(C(=O)C=C2)=O


Formula $C_{21}H_{18}F_5NO_6$ M.Wt 475.36

Solubility Soluble in DMSO Storage Store at -20°C

General tips For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution: ship with blue ice. All other available size: ship with RT, or blue ice upon request.

Structure



Background

Ald-Ph-amido-PEG2-C2-Pfp ester is a noncleavable 2-unit PEG linker for antibody-drug-conjugation (ADC).

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA