

Product Data Sheet

Product Name: BDW-OH

Cat. No.: GC91659

Chemical Properties

Cas. No. 573976-22-4


Formula $C_{11}H_{10}N_4O_2S_2$ M.Wt 294.4

Solubility DMSO: Sparingly Soluble: 1-10 mg/ml Storage -20°C

General tips For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution: ship with blue ice All other available size: ship with RT, or blue ice upon request.

Structure



Background

BDW-OH is an active metabolite of the STING agonist and prodrug BDW568. [1],[2]

References:

[1].Zhao, J., Tang, Z., Selvaraju, M., et al. Cellular target deconvolution of small molecules using a selection-based genetic screening platform ACS Cent. Sci. 8(10)1424-1434(2022).

[2].Tang, Z., Zhao, J., Li, Y., et al. Structural evaluations of a selective human STINGA230 agonist and its use in macrophage immunotherapies ACS Med. Chem. Lett. 15(5)653-658(2024).

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA