

Product Data Sheet

Product Name: KDdiA-PC

Cat. No.: GC43997

Chemical Properties

Cas. No. 439904-34-4

SMILES CCCCCCCCCCCCCCCC(OC[C@@H](OC(CCCCCCCC/C=C/C(O)=O)=O)=O)COP(OC[CN+](C)(C)C)([O-])=O


Formula C₃₆H₆₆NO₁₁P M.Wt 719.9

Solubility DMF: .5 mg/ml, DMSO: .5 mg/ml, Ethanol: 30 mg/ml, PBS (pH 7.2): 5 mg/ml Storage Store at -20°C

General tips For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution: ship with blue ice All other available size: ship with RT, or blue ice upon request.

Structure



Background

Oxidized low-density lipoprotein (oxLDL) particles contain low molecular weight species which are cytotoxic and pro-atherogenic. Many of these substances were recently isolated and purified from oxLDL, and identified as phosphatidylcholine species containing a fragmented, oxidized short-chain fatty acid remnant at the sn-2 position. KDdiA-PC is one of the most potent CD36 ligands among the oxLDL species. KDdiA-PC confers CD-36 scavenger receptor binding affinity to LDL at a frequency of only two to three KDdiA-PC molecules/LDL particle, and may be one of the more important structural determinants of oxLDL.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA

Product Data Sheet

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA