

Product Data Sheet

Product Name: Vedastotug

Cat. No.: GC79278

Chemical Properties

Cas. No. 2878367-20-3

Formula M.Wt 145.49

Solubility Storage Store at -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure

Vedastotug

Background

Vedastotug is a human monoclonal antibody targeting TNFRSF4/OX40/CD134, with its corresponding isotype control being Human IgG1 kappa, Isotype Control [1].

References:

[1]. International Nonproprietary Names for Pharmaceutical Substances (INN). WHO Drug Information, Vol. 39, No. 4, 2025

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA