

Product Data Sheet

Product Name: Dipotassium hydrogen phosphate (trihydrate)

Cat. No.: GB58250

Chemical Properties

Cas. No. 16788-57-1


Formula H₇K₂O₇P M.Wt 228.22

Solubility H₂O : 50 mg/mL (219.09 mM; Need ultrasonic) Storage Store at room temperature

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure



Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA