

Product Data Sheet

Product Name: Fmoc-Pro-Wang resin

Cat. No.: GA10050

Chemical Properties

Cas. No. N/A

Chemical Name 9H-fluoren-9-ylmethyl (2S)-2-formylpyrrolidine-1-carboxylate

SMILES C1CC(N(C1)C(=O)OCC2C3=CC=CC=C3C4=CC=CC=C24)C=O


Formula M.Wt

Solubility Storage Store at RT, Cool dry place

General tips For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution: ship with blue ice. All other available size: ship with RT, or blue ice upon request.

Structure



Background

Fmoc-Pro-Wang resin

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA