

Product Data Sheet

Product Name: Toddacoumaquinone

Cat. No.: GF01018

Chemical Properties

Cas. No. 142878-03-3

Formula C₂₃H₁₈O₇

M.Wt 406.38

Solubility

Storage 2-8°C, sealed storage, away from moisture, protect from light

General tips For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution: ship with blue ice. All other available size: ship with RT, or blue ice upon request.

Structure



Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA