

Product Data Sheet

Product Name: 6-Methoxy-1,2,3,4-tetrahydroquinoline

Cat. No.: GD17223

Chemical Properties

Cas. No. 120-15-0


Formula $C_{10}H_{13}NO$ M.Wt 163.22

Solubility Storage $-20^{\circ}C$

General tips For obtaining a higher solubility, please warm the tube at $37^{\circ}C$ and shake it in the ultrasonic bath for a while. Stock solution can be stored below $-20^{\circ}C$ for several months.

Shipping Condition Evaluation sample solution: ship with blue ice. All other available size: ship with RT, or blue ice upon request.

Structure



Background

Thalline (Thallin) is a tetrahydroquinoline derivative that can be used in chemical synthesis [1].

References:

[1]. Cromwell N H, et al. α , β -Diamino Ketones. II. 1 Reactions of Thalline and Open Chain Secondary Amines with α -Bromo- β -aminoketones[J]. Journal of the American Chemical Society, 1944, 66(3): 401-403.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA