

Product Data Sheet

Product Name: TN-783
Cat. No.: GC79411

Chemical Properties

Cas. No.


Formula $C_{16}H_{11}F_4IN_4O_2$ M.Wt 494.18

Solubility Storage Store at -20°C

General tips For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution: ship with blue ice. All other available size: ship with RT, or blue ice upon request.

Structure



Background

TN-783 is an orally active, brain-penetrant NLRP3 ($IC_{50} = 19.3$ nM) inhibitor. TN-783 enhances the effects of the GLP-1 receptor agonist semaglutide. TN-783 requires sustained targeted action, not permanent metabolic remodeling. TN-783 can be used for the study of obesity^[1].

In Vivo, TN-783 (0.5-50 mg/kg, p.o., twice daily, 14-66 days) can achieve safe (lean body mass preservation), effective, and reversible weight loss by reducing hypothalamic neuroinflammation^[1].

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA

Product Data Sheet

In Vitro, TN-783 (2 h) can effectively inhibit the secretion of IL-1 β (IC₅₀ = 45.9 nM) by THP-1 cells, under Lipopolysaccharides (LPS)/Nigericin stimulation^[1]. TN-783 (0.3 nM-6 μ M, 1 h) can dose-dependently inhibit the secretion of IL-1 β (IC₅₀ = 5.17 nM) and IL-18 (IC₅₀ = 12.3 nM) in THP-1 cells induced by Palmitic acid -BSA complex^[1].

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA