

Product Data Sheet

Product Name: TLR7/8 agonist 4

Cat. No.: GC65915

Chemical Properties

Cas. No. 2388520-33-8

Formula $C_{18}H_{24}N_6$

M.Wt 324.42


Solubility

Storage Store at -20°C

General tips For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT, or blue ice upon request.

Structure



Background

TLR7/8 agonist 4 (compound 41) is a potent **TLR7/8** agonist. TLR7/8 agonist 4 has anti-cancer activity^[1].

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA