

Product Data Sheet

Product Name: H-Arg-pNA.2HCl

Cat. No.: GA11100

Chemical Properties

Cas. No. 40127-11-5

Chemical Name (2S)-2-amino-5-(diaminomethylideneamino)-N-(4-nitrophenyl)pentanamide; dihydrochloride

SMILES C1=CC(=CC=C1NC(=O)C(CCCN=C(N)N)N)[N+](=O)[O-].Cl.Cl

Formula $C_{12}H_{18}N_6O_3 \cdot 2HCl$

M.Wt 367.3


Solubility

Storage

General tips For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution: ship with blue ice. All other available size: ship with RT, or blue ice upon request.

Structure



Background

H-Arg-pNA.2HCl

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA