

Product Data Sheet

Product Name: Dextran-CY5.5 (MW 10000)

Cat. No.: GC80373

Chemical Properties

Cas. No.

Formula M.Wt 10000

Solubility Storage Store at -20°C, protect from light, stored under nitrogen

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure

Dextran-CY5.5 (MW 10000)

Background

Dextran -CY5.5 (MW 10000) is a CY5.5 labeled Dextran conjugate. CY5.5 is a Cyanine dye (The excitation wavelength is usually 673 nm and the emission wavelength is 707 nm) that can be used for the labeling of proteins, antibodies and small molecular compounds. Dextran has an inhibitory effect

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA

Product Data Sheet

on thrombocyte aggregation and coagulation factors and is used as a plasma volume expander [1].

In Vitro, CY5.5-Dextran has a wide range of applications in biomedical research, such as studying vascular permeability and hemodynamics; labeling cells and then detecting cell migration and distribution using fluorescence microscopy or flow cytometry; serving as a carrier for studying drug distribution and metabolism in vivo; and for fluorescence imaging of tumor angiogenesis and metastasis.

References:

[1]. Vera DR, et al. Cy5.5-DTPA-galactosyl-dextran: a fluorescent probe for in vivo measurement of receptor biochemistry. Nucl Med Biol. 2005 Oct;32(7):687-93.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA