

Product Data Sheet

Product Name: Ald-Ph-amido-PEG4-C2-acid

Cat. No.: GC38346

Chemical Properties

Cas. No. 1309460-27-2

SMILES O=C(O)CCOCCOCCOCCOCCNC(C1=CC=C(C(=O)C=C1)=O


Formula $C_{19}H_{27}NO_8$ M.Wt 397.42

Solubility Soluble in DMSO Storage Store at -20°C

General tips For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution: ship with blue ice. All other available size: ship with RT, or blue ice upon request.

Structure



Background

Ald-Ph-amido-PEG4-C2-acid is a noncleavable linker used for the antibody-drug conjugate (ADC).

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA