

Product Data Sheet

Product Name: Anti-Mouse VISTA Antibody (8G8)

Cat. No.: GC81495

Chemical Properties

Cas. No.

Formula

M.Wt

Solubility

Storage Store at -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure



Background

Anti-Mouse VISTA Antibody (8G8) reacts with mouse VISTA. Anti-Mouse VISTA Antibody (8G8) is an agonistic antibody with immunosuppressive properties. Recommend Isotype Controls: Polyclonal

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA

Product Data Sheet

Armenian hamster IgG, Isotype Control [1].

References:

[1]. Mohamed A ElTanbouly, et al. VISTA is a checkpoint regulator for naïve T cell quiescence and peripheral tolerance. Science. 2020 Jan 17;367(6475):eaay0524.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA